Operating Instructions

A close-up photograph of the Ivoclar PrograScan PS7 laboratory scanner. The device is white with a green base and is illuminated with a green light. Two dental models are mounted on the scanner's platform. The Ivoclar logo is visible on the top of the scanner, and the model number 'PrograScan PS7' is printed on the front. The background is a blurred laboratory setting.

PrograScan[®]
PS7

Laboratory
scanner

PrograScan® PS7

[en]	Operating Instructions	4
	Laboratory Scanner	
[de]	Gebrauchsinformation	25
	Laborscanner	
[fr]	Mode d'emploi	45
	Scanner de laboratoire	
[it]	Istruzioni d'uso	65
	Scanner da laboratorio	
[es]	Instrucciones de uso	85
	Escáner de laboratorio	
[pt]	Instruções de Uso	105
	Scanner de Laboratório	

Table of contents

1	About this document	6
1.1	Signs and symbols	6
1.1.1	Warnings and additional information	6
1.2	Safety labels on the scanner	7
1.3	Marking	7
1.4	Customer support	7
2	Safety	8
2.1	Intended use	8
2.1.1.	Potential improper use	8
2.2.	Intended environment	9
2.3.	Operator's obligations	9
2.4.	Personnel qualifications	9
2.5.	Staff obligations	10
2.6.	Personal protective equipment	10
2.7.	Safety information for the individual operating phases	10
2.7.1.	Transportation and assembly	10
2.7.2.	Installation and operation	10
2.7.3.	Improper maintenance, modifications and repairs	11
3	Scanner overview and features	12
3.1	Scanner overview	12
3.2.	Connection panel	13
3.3.	Power supply unit	13
3.4.	USB cables	13
3.5.	Accessories	14
4	Function and assembly	16
4.1.	Unpacking the unit	16
4.1.1.	Scope of delivery	16
4.2.	Choosing the installation site	17
4.3.	Lifting points	17
4.4.	Place of installation	18
5	Installation	19
5.1.	Installing the scanner	19
5.2.	Installing the computer	19
5.3.	Installing the software	19
6	Operation	20
6.1.	ID recognition	20
6.2.	Optical Status Display	20
6.3.	Switch the unit off	20
7	FAQ – troubleshooting	21
8	Cleaning and maintenance work for the user	22
8.1.	Recalibration	22
8.2	Cleaning	22
9	Disposal	23
10	Technical data	24
10.1	Product drawings	24
11	Declaration of conformity	125

1 About this document

WARNING! Improper use of the product.

Risk to the patient and the user.

These Operating Instructions contain important information on the safety of the product: Please read the Operating Instructions thoroughly prior to installation and operation.

This optical 3D scanner (henceforth "scanner") is a product which guarantees the highest level of operational safety with state of the art technology. However, the scanner can be dangerous if used improperly, by untrained personnel or in a manner other than intended. Possible dangers are detailed in chapter "About this document" as well as being indicated in safety instructions throughout these Operating Instructions.

These Operating Instructions explain how to safely work on and with this scanner. The safety instructions it contains must be observed. These Operating Instructions contain the relevant information for the assembly, initial start-up, operation, maintenance and disposal of the scanner. The purpose of these Operating Instructions is to ensure risk-free operation. Observing the instructions in these Operating Instructions helps to avoid dangers and prevent damage to the product.

All persons must have these Operating Instructions at hand when working on and with the scanner and must observe the relevant information and notes. These Operating Instructions must always be complete and clearly legible. Ivoclar Vivadent AG is not liable for technical or typographical errors in these Operating Instructions, nor is any liability assumed for damages that are directly or indirectly attributable to the delivery, performance, or use of these Operating Instructions.

1.1 Signs and symbols

1.1.1 Warnings and additional information

Warnings

Warnings are used in these Operating Instructions to warn you of a risk of injury to persons and/or damage to property.

1. Please always read and observe these warnings.
2. Follow all measures marked with the warning symbol and warning word.

Depending on the severity and the probability of the danger, the following warning levels are distinguished:

Warning symbol	Warning word	Danger level	Consequences resulting from non-observance
	Danger	Immediate imminent danger	Death, severe personal injury
	Warning	Potential danger	Death, severe personal injury
	Caution	Potential danger	Slight personal injury
	Notice	Potential danger	Damage to property

Explanation of the structure of a warning:

WARNING! Name of source of danger, cause of hazard or type of risk.

- Consequences of non-compliance with instructions.
- Action in order to avoid danger.

Example of a warning:

DANGER! Residual energies may still be stored in the scanner, which may be unexpectedly released during maintenance.

Death or severe personal injury due to electric shock.

- Ensure that work is carried out only by qualified personnel and avoid any method of operation that could impair the safety of the scanner.
- Ensure that unauthorized persons cannot perform any work on the scanner.

Additional information

Symbol	Meaning
	Additional information, e.g. for better understanding, for simplifying workflows or for further information.

1.2 Safety labels on the scanner

Safety labels indicating possible dangers/residual dangers are mounted on the scanner.

The safety labels on the scanner must be followed under all circumstances. If they become faded or damaged during the scanner's lifetime, they must be replaced immediately.

Their readability and completeness must be checked at regular intervals. The moment the labels are not immediately recognizable and understandable at first glance, the optical 3D scanner must be taken out of operation until new labels have been attached.

The pictograms for warnings, prohibitions and commands on the scanner and their meanings:

Symbol	Description
	Danger of crushing This pictogram is located in the area of the swivel arm and is designed to warn of a possible crushing hazard between scanner and swivel arm.

1.3 Marking

The scanner can be uniquely identified by means of its identification plate. CE mark in accordance with:

- EMC Directive 2014/30/EU
- Low voltage directive 2014/35 EU
- RoHS Directive 2011/65/EU

1.4 Customer support

In case of technical problems, please contact your customer support or Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com.

2 Safety

These Operating Instructions are structured according to the applicable EU regulations and contain safety instructions. The company operating the scanner is responsible for ensuring that operating personnel receive the necessary safety information and that they have read and understood the Operating Instructions before using the scanner. It is the responsibility of the individual to comply with the safety instructions. The warnings used in these Operating Instructions do not include all safety instructions that must be followed. They only supplement general safety instructions and procedures. Ivoclar Vivadent AG accepts no liability for any consequences caused by the violation of general safety requirements or by the violation of design, production and usage safety standards. This chapter contains a general introduction to the safety instructions. It also contains important information on accident prevention.

Even if the greatest care is taken when designing and building the scanner and if all safety aspects are taken into account, residual dangers may still exist, which have been evaluated by means of a risk assessment. All residual risks and warnings from the risk assessment are listed in this chapter. The scanner complies with the latest technology and with the applicable health and safety regulations. Any statement, information or recommendation in this document does not constitute a promise of characteristics of any kind, either express or implied.

Nevertheless, the following dangers can occur in case of incorrect operation or misuse:

- Risk of fatal or personal injury to operator or third party
- Damage to the scanner or material assets of the operator
- Risk to the efficiency of the scanner

2.1 Intended use

The intended use of the scanner covers the three-dimensional, optical measurement of human jaw models in combination with the supplied software.

Only impressions and plaster casts of jaws as well as articulators with such casts may be used (1500 g maximum). Dimensions W x L x H: 200 mm x 200 mm x 150 mm. Only accessories specified by Ivoclar may be used. The scanner is designed to be used in dental laboratories and mounted at an ergonomic height (observe the health and safety regulations in your country).

Any other use requires the written consent of Ivoclar Vivadent AG. Incorrect use can lead to danger to persons and damage to the scanner. The scanner must be loaded manually, it does not function independently. For this reason the product may only be operated after the operator has read and understood these Operating Instructions. Furthermore, liability and warranty claims are invalid if the scanner is used in a manner other than intended. The scanner may only be operated in accordance with the operating conditions specified in these Operating Instructions.

2.1.1 Potential improper use

The information contained in these Operating Instructions describes the product's properties without guaranteeing them. No liability will be assumed for damage caused by:

- Improper use of the scanner
- Ignoring these Operating Instructions
- Unauthorized changes to the scanner
- Working on and with the scanner incorrectly
- Repairs which are carried out incorrectly
- Unauthorized, improper modification of operating parameters
- Adverse events resulting from foreign objects and force majeure

2.2 Intended environment

The conditions for the scanner's operation are limited as follows:

- Use only indoors
- Ambient temperature: 5 °C to 35 °C
- Relative humidity:
 - At temperatures up to 31 °C: 20% – 80%
 - At temperatures up to 35 °C: Maximum 20% – 50%, gradually decreasing
 - Not suitable for humid areas
- No flickering light (e.g. fluorescent tubes; moving blinds; shadows from trees in the wind)
- Maximum ambient brightness 500 lux; no direct sunlight in the scanning area
- Clean and dust free environment
- Install horizontally on a stable, vibration-absorbing workstation
- Maximum operating altitude: 3000 m above sea level

2.3 Operator's obligations

The operator is responsible for the safe operation of the device.

- Ensure that only instructed, trained or qualified persons work with the scanner.
- Ensure that every person working with the scanner for the first time has read and understood these Operating Instructions.
- Pay special attention to the safety instructions in these Operating Instructions and on the scanner.
- Keep these Operating Instructions in the scanner area and pass it on to new personnel.
- Make sure that nobody works on the scanner who does not have sufficient technical knowledge.
- Ensure that the operating personnel is instructed at regular intervals about safe scanner operating practices.
- Consider the health and safety regulations in your country.
- Ensure that during service, maintenance or cleaning work, safeguards are in place to prevent the scanner from being started up without authorization.
- Only use original replacement parts approved by the manufacturer. Otherwise guarantee regulations and product liability risks could be affected.
- Repairs may only be carried out by a certified service provider.
- Ensure that the labels attached to the scanner when delivered remain clearly legible. Missing or damaged warning signs must be replaced immediately.

2.4 Personnel qualifications

These Operating Instructions are intended for the personnel who operate and maintain the scanner. The operating and maintenance personnel must be selected by the operator and must fulfil the following requirements:

- Basic technical skills
- Read and understand these operating and maintenance instructions.
- The person must be strong enough to lift the scanner and move it to its final location.
- Operatives must not be visually impaired.

The following measures must be carried out by the operator to ensure appropriate knowledge of the scanner's operation:

- Product training
- Regular safety instruction
- Inspection, maintenance, cleaning and repair may only be carried out by technical specialists with product-specific training

2.5 Staff obligations

- Observe the safety instructions and operating instructions contained in these Operating Instructions at all times.
- Ensure that the scanner is only operated by suitably qualified personnel. These persons must be familiar with all safety instructions and with the corresponding measures relating to operation, maintenance and cleaning contained in these Operating Instructions and on the scanner.
- Before carrying out repair or maintenance work, disconnect the electrical and pneumatic supply to the scanner safely and properly safeguard it against being switched on.
- Disconnect the optical 3D scanner from the mains before cleaning, maintenance and repair work.

2.6 Personal protective equipment

The optical 3D scanner is designed such that persons working on the scanner do not require any additional protective equipment beyond the operator's required standard equipment (e.g. work clothing, work shoes, etc.).

- When carrying out cleaning and maintenance work, comply with the prescribed measures for the corresponding cleaning agents (e.g. gloves for cleaning agents, splash guard, etc.).

2.7 Safety information for the individual operating phases

2.7.1 Transportation and assembly

- Only transport the scanner if it is properly packed and secured.
- When handling the scanner, pay attention to the risk of bending the scanner. Especially when putting down the scanner after it has been transported.
- The scanner is designed to be mounted at an ergonomic height (check the health and safety regulations in your country). Ensure that the base underneath the scanner is capable of bearing at least twice the weight of the scanner.
- Pay attention to trip hazards when cables are laid on the floor. Always lay the cables along walls. The cables must not be tensioned under any circumstances or left dangling.
- Choose the scanner's location such that possible vibrations during the scanning operation can be absorbed without any problems, i.e. the load capacity of the base (workbench, table, etc.) must be at least twice the value specified in the technical specifications.
- When setting up the product, pay attention to the ambient light. Any change of the ambient light e.g. due to shadows from trees or blinds must be avoided during the operation of the scanner. In addition, avoid fluorescent tubes or direct sunlight.
- Do not use the scanner in damp rooms.
- Ensure that no liquids come into contact with the scanner. Do not place any vessels containing liquids in the vicinity of the scanner.

2.7.2 Installation and operation

- The operating and maintenance personnel must be instructed regularly on how to operate the system as well as on its dangers. It is mandatory to read the Operating Instructions.
- Ensure that only qualified personnel installs the system.
- Every time you start up the scanner, ensure that no persons are in the pivoting range of the swivel arm since it performs a reference run when the scanner is switched on.
- Every time you start up the scanner, ensure that no persons are in the danger zone. In case of danger, the scanner must be disconnected from the mains immediately.
- Keep children away from the scanner. Only qualified personnel may work with the scanner.
- Only operate the scanner in perfect and clean condition.
- Ensure that power supply cables and feeder cables are not damaged or crushed by other objects. Report faults or damage to the operator immediately.
- Do not remove, modify, bridge or bypass any protective, safety or monitoring device.

- The swivel arm moves during operation. Do not reach into the danger area as this may result in injury to your fingers.
- When operating the scanner inserts, pay attention to the possible danger of crushing.
- If the scanner is not going to be used for an extended period of time, disconnect it from the power supply.
- Do not place objects or tools on the scanner.
- Keep the area around the scanner clean and tidy. Remove all objects and containers that are not required for the work process from the scanner area.
- Do not climb the scanner.

2.7.3 Improper maintenance, modifications and repairs

Retrofitting additional devices from third-party manufacturers is not permitted, nor is the modification of protective devices without consultation with Ivoclar Vivadent AG or the manufacturer of the respective additional device.


- Ensure that only personnel of the Ivoclar Vivadent AG customer service repairs the system.
- Opening the scanner is strictly prohibited. Ensure that the scanner is only opened by Ivoclar service personnel.
- Only clean the scanner with a dry cloth.
- Do not convert or modify the scanner without written approval from Ivoclar Vivadent AG.
- Disconnect the scanner from the power supply when carrying out repairs.
- Use only original spare parts.
- After maintenance or repair work, ensure that all protective devices/casings are mounted before the scanner is started.

3 Scanner overview and features

This optical 3D scanner is designed for use in the dental field to measure and digitize human jaw models in three dimensions. Using the "structured-light scanning" principle, a stripe pattern of light is projected onto an optical mirror. The deflection of the light pattern caused by the structure of the model is registered and processed by the four high-resolution cameras of the scanner unit, and subsequently digitized into a virtual, three-dimensional model.


3.1 Scanner overview

The scanner's main components are as follows:

3.2 Connection panel

The connection panel is located on the rear of the scanner and is laid out as follows:

3.3 Power supply unit

NOTICE! Only use with the supplied power supply unit and power cable!

- Otherwise guarantee regulations and product liability risks could be affected.
- Only use with the supplied power supply unit and power cable.
- Only connect the power supply unit to a socket with protected contacts.

Technical data of the power supply unit:

Manufacturer	MEAN WELL
Type	GST90A 24 – P1M
Input	90 ~ 264VAC 127 ~ 370VDC 47 ~ 63Hz
Power	90 W
Output	24 V, 3.75 A

3.4 USB cables

NOTICE! Check the power of the USB cables.

- Using two USB cables could result in lower performance.
- Connect the USB cables to two separate USB 3.0 ports. The full bandwidth is required for each USB cable.

Use the supplied USB cable. Should you use your own USB cable, make sure the following requirements are complied with:

Shielded, full-featured 5 Gbit USB-C cable, maximum 1.4 m length. In the case of extensions over 1.4 m, use an active USB extension.

3.5 Accessories


NOTICE! Do not use non-original parts

Otherwise guarantee regulations and product liability risks could be affected.

– Use original parts only.

The scanner is delivered with standard equipment. Accessories for certain functions may be purchased.

Standard equipment	Item number	Description
	760800	Calibration body autodetection
	760801	Articulator base autodetection
	760802	Multi-die holder (9) autodetection
	760804	Model mounting Blu Tack autodetection
	760885	Impression holder autodetection
Accessories	Article numbers	Description
	760803	Model mounting screw system auto detection
	760812	Model fixator rod
	760805	3D impression holder autodetection

	760806	Multi-die holder (14) autodetection
	760807	Quarter tray holder autodetection
	760808	Adesso SplitCast holder autodetection
	760809	SAM/Axio SplitCast holder autodetection
	760810	Protar SplitCast holder autodetection
	760811	Quicksplit SplitCast autodetection
	760813	Articulator calibration object

4 Function and assembly

NOTICE! Transportation in inadequate or faulty packaging.

Damage to the scanner.

- Ensure that the scanner is transported exclusively in the packaging provided.
- Do not expose the scanner to moisture or heat at any point during transport.
- When stored for extended periods, keep the scanner in a dry place in its original packaging. This should prevent corrosion and contamination.

Make sure that the reference points of the calibration plate are not damaged as otherwise calibration is no longer possible.

- Check if the packaging is undamaged. If this is not the case, document the damage and contact your dealer.
- Pay attention to the warnings on the packaging when transporting the scanner.
- Transport or store the scanner well protected, dry and in upright position in its original packaging.
- Ensure that the storage provisions in the specifications are complied with.

4.1 Unpacking the unit

Dispose of the packaging in an environmentally friendly way (paper into paper recycling, plastic into plastic recycling, etc.). We recommend keeping the packaging in case of returns.

- Open the top of the packaging and remove the foam cover with the accessories.
- Pull the scanner out of the packaging via the two notches on the sides.
- Remove the holder of the swivel arm. Inside the holder are the accessories delivered with the scanner.
- Check that all parts included in the delivery form are present.

4.1.1 Scope of delivery

- 1 x PrograScan PS7
- 1 x Calibration body autodetection
- 1 x Articulator base autodetection
- 1 x Multi-die holder (9) autodetection
- 2 x Model mounting Blu Tack autodetection
- 1 x Impression holder autodetection

4.2 Choosing the installation site

WARNING! Scanner exposed to direct sunlight or strong extraneous light.

Fire hazards and damage to the scanner due to overheating and unwanted reflections on the scanner's monitor and sensor technology.

- Do not expose the scanner to direct sunlight and avoid direct proximity to heat sources (e.g. radiators, other electrical appliances, fireplaces, etc.).
- Ensure generous air circulation around the scanner.

Since the environment has a significant influence on the scanner, you should be careful when selecting an installation site and observe the following points:

- Avoid rapidly changing lighting conditions such as moving shadows and flickering lighting (fluorescent LED).
- Do not expose the scanner to direct sunlight.
- Avoid direct proximity to heat sources (e.g. radiators, other electrical appliances, fireplaces, etc.).
- Ensure that the installation site is level and free from vibrations and contamination.
- Ensure that no extraordinary exposure to dust or toxic or corrosive gases and vapours or to impermissible heat effects occurs at the installation site.
- Select a mounting or installation site where liquids or objects cannot enter the scanner under any circumstances (e.g. condensation, roof leaks, spills, etc.).
- Ensure that the installation site is equipped with the necessary power and network connections.
- Only use and install the scanner indoor.
- Ensure that the installation site is free of increased electromagnetic radiation.
- Select a mounting or installation site which prevents children from touching the appliance and its connections unattended.
- The installation site must allow for the safe positioning of all connected cables.
- Ensure that power supply cables and feeder cables are not damaged or crushed by other objects.
- Ensure that the mains plug of the power supply unit is always accessible.

4.3 Lifting points

The two lifting points (1) and (2) are for carrying the scanner. Hold the scanner at these two points with your left and right hand. The scanner weighs 11–13kg and can be carried by one person, with the open side of the scanner facing away from the body.

4.4 Place of installation

WARNING! Scanner exposed to direct sunlight or strong extraneous light.


Fire hazards and damage to the scanner due to overheating and unwanted reflections on the scanner's monitor and sensor technology.

- Do not expose the scanner to direct sunlight and avoid direct proximity to heat sources (e.g. radiators, other electrical appliances, fireplaces, etc.).
- Ensure generous air circulation around the scanner.

Make sure to run cables along walls. The power pack must be connected to a socket with a protected contact.

The following diagram is a schematic representation only which serves as a visualization of the overall system. Pay attention to an ergonomic arrangement as well as to the drawn-in minimum distances of the product (grey).

5 Installation

NOTICE! Only use with the supplied power supply unit and power cable!

- Otherwise guarantee regulations and product liability risks could be affected.
- Only use with the supplied power supply unit and power cable.
- Only connect the power supply unit to a socket with protected contacts.

5.1 Installing the scanner

NOTICE! Increased electromagnetic radiation.

- Damage to the scanner.
- Ensure that the installation site is free of increased electromagnetic radiation.

NOTICE! Check the power of the USB cable.

- Using two USB cables could result in lower performance.
- Connect the USB cables to two separate USB ports. The full bandwidth is required for each USB cable.

After setting up the scanner, follow these instructions on how to install it:

- Connect the scanner and PC to the USB ports using the supplied USB cables. If the USB cables are too short, they can be extended with special active USB extensions. For appropriate products, please contact Ivoclar Vivadent AG.
- Connect the power supply unit to the scanner and the mains using the power cable.

5.2 Installing the computer

- The computer's instruction manual is pre-installed on the operating system.
- Check the computer's minimum technical requirements (Windows 10 or higher) to ensure the proper functioning of the scanner.
- Make sure that all drivers are up to date.
- Install the computer according to the computer manufacturer's installation instructions.

5.3 Installing the software

- For quick assistance, please have your access data ready before contacting our support. Our support works with TeamViewer (www.teamviewer.com).
- Check the computer's minimum technical requirements (Windows 10 or higher) to ensure the proper functioning of the scanner.
- Make sure that all drivers are up to date.
- Install the computer according to the computer manufacturer's installation instructions.

6 Operation

CAUTION! The swivel arm moves when the scanner is switched on.

Danger of crushing.

- Do not reach into the danger area when switching on the scanner.

Make sure that the reference points of the calibration plate are not damaged as otherwise calibration is no longer possible.

Prerequisite:

- The scanner has already been set up and the installation has been carried out.
- The scanner is switched on.

Follow the steps below:

- Start the scan software by selecting the following icon:
- Place the calibration plate in the scanner.
- Select scanner calibration from the Options menu and start the calibration.
- As shown in the menu, you can create a new case after successful calibration.
- To scan an object, place it in the scanner and press the scan menu.

For more details, please download the PrograScan software manual. It is available in the software.

6.1 ID recognition

Each accessory has its own ID, which enables the scanner to determine when it is in the unit. This allows the corresponding functions to be carried out.

6.2 Optical Status Display

The LED display shows the basic operating statuses of the software.

- Green and blue in a slow colour gradient from left to right => Power on, Standby, Ready
- White light from left to right - Working process running
- Pulsating white light - No connection to the software
- Rainbow - Calibration in progress
- Progress bar - Scanning progress

6.3 Switch the unit off

To switch off the product, press the On/Off button for at least 1 second.

7 FAQ – troubleshooting

Scanner does not start	<ul style="list-style-type: none">– Check the power supply if the scanner does not light up at start-up.– Contact our customer service. See chapter "Customer support".
The scanner is not recognized by the PC.	<ul style="list-style-type: none">– Make sure you use high-quality USB cables with a maximum length of 1.4 m (or an active USB extension from CADstar). Test by changing the cable.– Check the required USB 3.0 specification of the USB port.– Check that all drivers are up to date.– Check that the LED on the scanner is lit.– Restart the PC and the scanner– Contact our customer service. See chapter "Customer support".
The scan result shows streaks	Pay attention to light changes during the scan. Light changes can lead to streaks on the scan result. This can be caused by blinds moving during the scan or sunlight and flickering light. You can test the lighting conditions with the help of a slow-motion recording from your smart phone.
The scan result shows holes	<ul style="list-style-type: none">– Perform an integration using the integration tool.– If there are dark spots on the object to be scanned, use a scan spray.

8 Cleaning and maintenance work for the user

To ensure long product life and scan accuracy, it is necessary to regularly clean and calibrate the product.

8.1 Recalibration

- **The calibration plate must not be cleaned. Touching the reference points located on the calibration plate can cause the calibration accuracy to decrease.**
- **Make sure that the reference points of the calibration plate are not damaged as otherwise calibration is no longer possible.**
- **In case of a malfunction that cannot be resolved by restarting the scanner or workstation, contact your customer support.**

To guarantee consistent scan results, we recommend you calibrate the scanner on a weekly basis using the supplied calibration plate.

The scanner should also be calibrated after transportation or if the scanner is moved. Instructions for recalibrating your scanner can be found in our PrograScan software manual. Initiate the scan process using the software and follow the instructions step by step. If you are unsure, please contact our customer support.

8.2 Cleaning

DANGER! Accidentally switching on the scanner during cleaning.

Death or severe personal injury due to electric shock.

- Before cleaning, take the optical 3D scanner out of operation and safeguard it from being started without authorization.
- Turn the scanner off using the switch on the rear panel to eliminate any risk of injury caused by accidentally switching on the scanner.

NOTICE! Scanner exposed to water, steam or cleaning agent.

Damage to the electrical parts of the scanner.

- Ensure that no water, steam or cleaning agents penetrate the electrical scanner parts during cleaning.

NOTICE! Improper cleaning of the optical sensor. The optical sensor located in the upper part of the scanner is extremely sensitive.

Damage to the optical sensor.

- Do not clean the optical sensor system located in the upper part of the scanner.
 - Do not use harsh cleaning agents, as they can also damage the scanner.
-
- Ensure that the scanner is only cleaned by people who are trained and who have been informed of the scanner's potential dangers.
 - Stop the scanner and disconnect it from the power supply for the duration of the cleaning.
 - Vacuum or wipe the scanner with a dry cloth.
 - Pay attention to the optical components when cleaning. These must not be adjusted.
 - To prevent damage to the optical sensors located in the upper part of the scanner, they may only be dusted with a microfibre cloth.

9 Disposal

- In case of disposal, the scanner can be returned to the dealer or manufacturer.
- Please note that the scanner is an electronic device that is exclusively for use in commercial or industrial areas.
- Disposal using public waste management companies is therefore not possible.
- Please contact our customer support for disposal.

10 Technical data

General description	
Manufacturer	Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan, Liechtenstein
Registration	2023
Weight	13.6 kg
Dimensions (W x H x D)	396.0 x 436.5 x 380.5 mm
Environment	
Operating temperature range	18 °C to 35 °C
Storage temperature range	-15 °C to 50 °C
Altitude	0 to 3000 m
Pollution degree	2
Electrics	
Supply voltage	100–240 VAC
Mains frequency	50–60 Hz
Performance power supply unit	90 Wh
Voltage scanner	24V DC +-10%
Power consumption operation	45 Wh
Power consumption stand-by	1 Wh
Camera	2x (4x 8 Mpix)
Type of USB connection	USB 3.0 TYP C (5 Gbit/s)
Number of USB ports	2
Scanner	
Dimensions of measurement area (W x H x D)	2x 100 x 100 x 100 mm
Maximum object size	300 x 200 x 150 mm
Maximum object weight	2 kg
Absolute precision	<5 μ
Repeat precision	<2 μ
Output format	STL, PLY

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11 Declaration of conformity

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Template ID: TEFO-04701-EN Version:1.0, Valid as of: 14 Sep 2022 11:23:21 (GMT+02:00)
 Template title: 02-602 Declaration of Conformity (non MD)
 Module: Device Registration and Facility Registration (Non-EU)
 Module owner: Head of Department International Registration

We herewith declare under our exclusive responsibility that the below mentioned products meet the provisions of the following EC Council Directives and its implementation in national law. All supporting documentation is retained on the premises of the manufacturer and, where applicable, the notified body.

Product	PrograScan PS7-G2
Document-ID	LL 11962946
Document Version	1.0

Legal manufacturer	
	Ivoclar Vivadent AG Bendererstrasse 2 9494 Schaan/Liechtenstein www.ivoclar.com

Declaration of Conformity Information	
Directives for equipment	<input checked="" type="checkbox"/> 2014/35/EU Low voltage <input checked="" type="checkbox"/> 2014/30/EU Electromagnetic compatibility <input checked="" type="checkbox"/> 2011/65/EU RoHS
Standards for equipment	<input checked="" type="checkbox"/> EN IEC 62368-1:2020 + A11:2020 <input checked="" type="checkbox"/> EN 61326-1:2013 (Group 1, Class B) <input checked="" type="checkbox"/> EN IEC 63000:2018
Notified Body Address	n/a
Place and date of issue	Schaan, 2023-06-22
Valid until	2028-06-21

Document Control		
Name	Date	Signature
Approver (Director) Gottfried Rohner	22.06.2023	
Approver (CTO): Dr. Thomas Hirt	22.06.2023	

Date information prepared: 2023-06 / Rev.0

ivoclar.com