ORBIS

ORBIS Compomer

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Instructions for use

In accordance with DIN EN ISO 4049

ORBIS Compomer is a universal light-curing radiopaque compomer restorative material for filling cavities in both the anterior and posterior areas. It allows a simplified and rapid filling technique as well as combines the fluoride release of glass ionomers with the physical durability and high aesthetics of composites.

ORBIS Compomer is used in combination with a suitable bonding agent.

Indications:

- filling deciduous teeth
- class V fillings (cervical caries, root erosions, V-shaped defects)
- facetting of discoloured anteriors
- class III fillings
- temporary class I and II fillings
- extended fissure sealing

Contraindications:

- If it is not possible to provide a dry gingival sulcus or the prescribed application technique, or there are known allergies to individual components (methacrylates, amines), the application should be avoided.
- ORBIS Compomer** is contraindicated for direct or indirect pulp capping as well as for occlusion-loaded, permanent fillings in secondary dentition.

Application:

Preparation/shade selection:

Clean the teeth to be treated with fluoride-free cleaning paste. To select the shade, compare the shade guide with original material (corresponds to the shade of an applied and thoroughly polymerised filling) with the cleaned tooth before anaesthesia (moisten shade guide and the tooth, select the shade in daylight). Secure a dry working field.

Cavity preparation:

In general, preparation should be minimal (adhesive filling technique) to conserve healthy tooth substance. Observe special preparation forms for extended fissure sealing and on deciduous teeth. Caries-free cervical defects do not require preparation. Clean the defect with either cleaning paste or pumice using a rubber cup or rotating brush.

Cleaning:

- Remove all residues from the cavity with the water jet.
- Avoid any contamination with blood or saliva after the cleaning.
- Remove excess water with a faint air jet, yet do not overdry the dentine. A moist dentine surface is necessary for wet-bonding.

Pulp protection:

A calcium hydroxide lining should be applied in close proximity to the pulp.

Application of a bonding:

The manufacturer’s instructions for use must be observed.

Application of ORBIS Compomer:

Apply and light-cure fillings more than 2 mm deep in layers. Polymerise each layer for 40 s with a conventional polymerisation device. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices.

Use transparent strips and crowns or light wedges. Ensure a tight marginal seal. Hold the tip of the light rod as close as possible to the surface of the filling. If the distance exceeds 5 mm the curing depth may be reduced to polymerisation in the range of the light cone only.

Incomplete polymerisation may lead to discoloration and pulpitis-like complaints.

Finishing and polishing of the filling can be conducted with cooling provided (e. g. fine or extra fine finishing diamonds, polishing discs) immediately after removing the shaping aids. As the last step, the margin of the filling or the entire tooth should be fluoridated.

Warnings, precautionary measures:

- ORBIS Compomer** contains Bis-GMA, diurethane-dimethacrylate, TEGDMA, BHT. **ORBIS Compomer** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Preparations containing eugenol or oil of cloves disturb the polymerisation of resin-based filling materials and should not be used in combination with resin-based filling materials.
- Contact with cationic mouth rinses as well as plaque indicators and chlorhexidine may lead to discoloration.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage:

Storage at 4 °C - 28 °C. Close syringes and bond bottles immediately after dispensing to avoid premature curing and/or evaporation of the solvent. Do not use after expiry date.

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**MP Medical Product GmbH**