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## GB

Caution: the product must be cleaned, disinfected and sterilised before first use as well as after any further use, according to the enclosed instructions.

### Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: HAGER & WERKEN GmbH & Co. KG  
Ackerstraße 1, 47269 Duisburg, Germany  
T +49 (203) 99269-0  
F +49 (203) 299283  
www.hagerwerken.de

#### Medical device

Product: Miratray 140°

REF 101 155 – REF 101 162

**Restriction for reprocessing:** max. 50 recycling

**Location of use:** none

**Storage and Transport:** No special requirements

Cleaning preparation: Remove impression tray of impression materials to remove surface contamination

**Cleaning:** Use a suitable cleaning agent according to the instructions of the manufacturer. After cleaning remove detergent thoroughly under running water. With cleaning device do not exceed 93 °C

**Disinfection:** Use suitable disinfectant solution according to the manufacturer's instructions. Remove after disinfection disinfectant thoroughly under running water. With disinfection device exceed 93 °C not

**Maintenance:** No special requirements

**Packaging:** Standardized packaging material for sterilization can be used (recommended)

**Sterilization:** Steam sterilization, 134 °C, holding time 5 min or Steam sterilization, 121 °C, holding time 15 min

**Inspection / functional test:** Visual inspection for damage, deterioration, deformation

**Storage:** No special requirements

The above statements have been validated by medical device manufacturers for the preparation of a medical device to its reuse as SUITABLE. The user is responsible for ensuring that the actual treatment performed with appropriate equipment, materials and personnel in the processing device achieves the desired results. This will normally require validation and routine monitoring are required. Likewise, any deviation from the instructions provided by the preparer carefully evaluated for their efficacy and possible adverse consequences.

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