ORBIS

LC Seal F

EN **DE** **DA** **SV** **NO** **FI**

PL **LT** **ET** **LV** **NL** **CS**

MD **CE** **0482** **REF** 309963

EN Instructions for use

Light cure pit and fissure sealant

Product description
ORBIS LC Seal F is a light cure single component and fluoride containing methacrylate based pit and fissure sealant; semi-opaque and white.

The white colour allows controlled application and facilitates visual observation, the semi-opaque colouration is more tooth like.

Etching Gel roughens the enamel surface and produces microfine retentions into ORBIS LC Seal F can flow. Thus, etching of the enamel increases adhesion of ORBIS LC Seal F in fissures, pits and on the tooth surface.

ORBIS LC Seal F meets the requirements of: DIN EN ISO 6874, class 2.

Indications/Intended use

Sealing of fissures and pits.

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

In singular cases, the material may cause a sensitizing reaction in patients with a hy-persensitivity to any of the ingredients. In these cases, the material should not be used.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhib-its the polymerization of the material. Neither store the material in proximity of eug-enol containing products, nor let the material allow coming into contact with materials containing eugenol.

Application

Isolation

Rubber dam is the recommended method of isolation (e.g. ORBIS rubber dam).

Enamel etching

The occlusal surfaces of the teeth to be sealed must be thoroughly cleaned. Apply Etch-ing Gel (e.g. ORBIS Etching Gel) in the area of the surface and keep the acid off from dentin. Leave the Etching Gel in place for 30 seconds and flush than thoroughly with water. Dry with air that is oil and water free. The etched enamel should have a chalky white appearance.

Etching procedure: It is essential that the etched enamel is not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

During the application of ORBIS LC Seal F, the areas should be kept dry with rubber dam.

Application of ORBIS LC SEAL F

1. Direct Application with Application Tips:

Screw off the cap of the syringe and substitute it by one of the supplied Application Tips.

Press cautiously ORBIS LC Seal F in a thin layer directly onto etched areas and – if necessary - spread with a suitable brush.

After use discard the Application Tip and attach the cap to the syringe again. Application Tips are for single use only due to hygienic reasons.

2. Application with a brush:

Place a small quantity of ORBIS LC Seal F onto the tooth and spread with a suitable brush.

ORBIS LC Seal F should be applied only to the etched areas. Remove excess of material.

By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm², cure each segment for **20 seconds**. After polymerization, adjust occlusion and inspect sealant for complete coverage. In case of incomplete coverage reapply ORBIS LC Seal F in the same manner.

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the material.
- Do not use any resin to adjust viscosity of the material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized material may have an irritant effect and may cause sensitization or allergic reaction against methacrylates in sensitive individuals.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Composition

Dimethacrylates, dental glass, photo initiators.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The European safety and clinical performance of the medical device can be found in the summary database on medical devices (EUADMED – https://ec.europa.eu/tools/euadmed).

Warranty

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, you exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental GmbH is not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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