
# TwinWave

USER MANUAL  
BENUTZERHANDBUCH  
MANUEL D'UTILISATION  
MANUAL DEL USUARIO  
MANUALE D'USO


## CONTENT

|                            |    |
|----------------------------|----|
| USER MANUAL.....           | 4  |
| BENUTZERHANDBUCH .....     | 24 |
| MANUEL D'UTILISATION ..... | 44 |
| MANUAL DEL USUARIO.....    | 64 |
| MANUALE D'USO .....        | 84 |

# 1. DEVICE CONTENTS


Tip  
Handpiece  
Battery pack  
Charging base  
Adaptor  
Rubber shield  
Protective fllasses

Tip

Handpiece

Battery pack

Charging base

Adaptor

Shield

Protection glasses

## 2. DEVICE DESCRIPTION

The device is a cordless, pen-style, LED light polymerisation device and fluorescence detector for use by dental professionals in dental offices or dental laboratories.

- Small size and lightweight ergonomic design.
- Compact cordless design and replaceable battery pack.
- Individually adjustable LED tip, can be rotated by 360°.
- LED tip design providing excellent intra-oral access.

### Indications

- For light-activated polymerisation of dental materials such as composites, luting cements and sealants using visible light.
- Detection of fluorescence of appropriate dental materials.

### Contraindications

The device is contraindicated for use in patients prone to photobiological reactions (including patients with solar urticaria or erythropoietic protoporphyria) or those currently undergoing treatment with photosensitising pharmaceuticals.

### Compatible materials

The device is designed to cure materials initiated with camphorquinone and/or other initiators absorbing violet light, 405–480 nm wavelength. See the polymer-based restorative material manufacturer's complete instructions for use for specific product compatibility and curing recommendations.

### 3. SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these instructions for use.

#### Safety

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Comply with all safety messages that follow this symbol in order to avoid possible injury.

#### Warnings

Never modify the device or any of its equipment. Any modifications may compromise safety and effectiveness.

#### Handpiece

- When using the fluorescence detector function, don't look into the light without protection; it's harmful to your eyes. In order to avoid the eyes coming into contact with the light directly, please assemble the protective glasses.
- Do not use a damaged device, e.g. if the glass cover of the LED tip is scratched, broken or missing.
- The eye protection shields will wear out over time. Replace the defective shield with a spare shield included in the package or available as an accessory.
- Only authorised technicians should repair the handpiece or battery pack.

#### Charging base

- Low voltages are present inside the charging base. Use only under dry conditions. Do not use if the charging base or handpiece is wet. Avoid a short circuit between the contact pads on the charging base. Only authorised technicians may repair the charging base.
- Do not use with voltages different from the range indicated on the charging base and power adaptor.
- Always make sure the handpiece, tips and battery pack are completely reprocessed and dry before inserting into the charging base.

#### Battery pack

- Prevent a battery short circuit during use and storage.
- Keep electrical contacts clean and dry.
- Do not remove the battery pack from the handpiece during operation.

## Application

- Do not use the device for intraoral illumination or dental transillumination. Excessive heat may develop, causing burns to mucosa or pulpal irritation.
- Insufficient material curing may lead to postoperative sensitivity and/or premature restoration failure.

## Precautions

This device is intended to be used only as specifically outlined in these instructions for use.

- Do not use a device which has not been properly reprocessed. Protect the device from gross contamination by applying the single-use barrier sleeve. Barrier sleeves are intended for single use only in order to prevent cross-contamination.
- Never aim the light directly at unprotected soft tissues, as this may cause injury or irritation. Do not aim the light at eyes. Light reflected from the tooth surface may also cause injury to eyes. Use the eye protection shields supplied with the device or protective glasses.
- Limit the action of the light to the area being treated.
- All dental curing lights cause a certain degree of heat development. Extended operation in areas near the pulp or soft tissues may result in severe damage. In these circumstances, do not cure for more than 10 seconds at a time without taking precautions such as air cooling.
- During extensive use (multiple curing cycles with 30 s or less between cycles), it is possible that the probe tip, which is an applied part, may reach a temperature of up to 50 °C. There should be no adverse effects resulting from short-term contact with intact skin or mucosa.
- Use only the power supply, power cord, charging base and battery supplied by the manufacturer. Use of any accessories other than those specified in the instructions for use may result in damage to the device and its components.
- Sterilising the device, components and accessories will cause component damage and may cause injury to the body. The eye protection shields may be autoclaved.
- Failure to follow recommendations for environmental operating conditions may result in injury to patients or users .
- Inspect the equipment before each use for worn, loose or damaged parts.
- The user may not service any parts except the O-ring attached to the coupling end of the light tips. Opening any of the components may result in unsafe operation and will void the warranty.
- This device must not be used in the presence of a flammable anesthetic gas mixed with air, oxygen or nitrous oxide.
- The user must not touch the patient and accessible charging base contacts simultaneously. Wear suitable protective eyewear, clothing, gloves and a mask. Protective eyewear is also recommended for patients.
- Do not spray disinfectant or other liquid directly onto the light, tips, battery, charging base, power supply or cord.

- Prevent liquids from entering the curing light body (handpiece), battery back and charging base.
- Do not place the system on or next to a radiator or other heat source. Excessive heat may damage the system's electronics.
- The tip, as well as the light emitted by the tip, must not come into contact with the patient's mouth.

#### Adverse reactions

- Prolonged unfiltered exposure to the light source may cause damage to the eye. (See Warnings).
- Prolonged contact with soft tissue may cause tissue injury or irritation. (See Warnings).
- Medical conditions such as solar urticaria, erythropoietic protoporphyria or cataract surgery may be aggravated by exposure to emitted light. (See Contraindications, Precautions).

## 4. STEP-BY-STEP INSTRUCTIONS


### Installation

1. Attach the tip to the handpiece by firmly pressing the tip into the handpiece opening while slightly rotating at the same time.
2. Adjust the LED tip: the LED tip can be rotated by 360°, thus the position of the light emission may be individually adjusted.
3. Insert the battery pack into the handpiece.

### Charging

1. Connect the charging base to the adaptor cord.
2. Make sure the power socket used for the power connector is accessible at all times in case an emergency disconnection is necessary.
3. Insert the battery into the charging base.


**Note:** The device comes with two batteries. It is recommended that the battery not being used be stored on the charging base so that it is fully charged when needed.

#### 4. Charging light

**Note:** No light when the battery is inserted indicates the battery is not in contact with the charging base.


- Solid orange: battery is charging
- Solid green: battery is fully charged

#### Operation at a glance

**Rotary switch:** The battery pack is also a rotary switch. You can choose different functions by operating it.

- Cure: curing light
- Check: fluorescence detector
- Off: power off

**Indicator lights:** Under the key, solid light or flashing indicates different functions.


- Solid green: curing light
- Slowly flashing green: fluorescence detector
- Quickly flashing green: battery has low power

**Note:** If the battery has low power, it must be replaced to continue to operate.

**Start key:** starts or disrupts working.

- Cure: press once for 10-second cycles, twice for 20-second cycles.
- Check: press once for 30-second cycles, twice for 60-second cycles.

**Note:** if there is no key operation for 3 minutes in standby mode, the device will automatically enter sleep mode in order to save power. Press the Start key to wake it up.


### Curing light radiometer

There is a light radiometer on the charging base to test the intensity of the curing light.


**Note:**

- The light radiometer can only be used to measure the curing light. For fluorescence detection light, the measured value is inaccurate.
- Before measuring, please turn the Rotary switch to Cure.
- Solid orange: light output is below 1000 mW/cm<sup>2</sup> and is inadequate (e.g. incorrect positioning, contaminated or scratched lens).
- Solid green: indicates irradiance of at least 1000 mW/cm<sup>2</sup>.


### Preparation before operation

1. Protect the LED tip from gross contamination by using a single-use barrier protection sleeve. Make sure that the lens is not blocked by crinkles or the sleeve's seam.

#### Cross-contamination

- Ensure the disposable polyethylene barrier protection sleeve has been correctly applied over the entirety of the LED tip.
  - The polyethylene barrier sleeve is not a substitute for cleaning and disinfection of the dental instruments.
2. Attach the eye protection shields supplied with the device in combination with barrier protection sleeves. Hold the tip over the opening of the shield and turn the tip against the shield by 90° into its final position.
    - Always make sure eye protection shields are securely fixed to the device to avoid accidental aspiration (press the eye protection shield firmly into the correct position). Always make sure the eye protection shields are properly mounted on the device without covering the light aperture.
    - Always make sure the eye protection shields are securely fixed to the tip to avoid accidental aspiration (press the eye protection shield firmly into the correct position).
    - Always make sure the light aperture is not covered by the eye protection shields.

### 3. Use suitable, light filtering safety glasses

#### **Eye damage caused by high-intensity light**

- Do not press the activation button until properly positioned intraorally.
- Ensure that everyone in the operating area (patients, operators, assistants) is wearing appropriate protective filtering eyewear.
- Do not look directly at the light while it is activated.

### Operation – curing light

1. Turn the rotary switch to Cure.
2. The LED tip should be positioned as closely as possible to the restoration. Avoid casting shadows by angulating the long axis accordingly. Set to the correct position.
3. Polymerisation: press the Start key to activate the light.
  - Press once for 10-second cycles.
  - Press twice for 20-second cycles.
4. To stop the curing light before the end of the cycle, press the Start key at any time.
5. Please refer to the instructions for use for the respective material and apply the curing times specified for 1000 mW/cm<sup>2</sup>.

### Operation – fluorescence detector

1. Turn the rotary switch to Check.
2. Check: Press the Start key to activate the light.
  - Press once for 30-second cycles.
  - Press twice for 60-second cycles.

## 5. CLEANING, DISINFECTION AND STERILISATION

### Cross-contamination

- Do not reuse single-use products. Dispose of them in accordance with local regulations.
- The barrier is designed for single use and must be disposed of after each use in accordance with local regulations. The barrier is not a substitute for cleaning, disinfection and sterilisation.

### Electrical short-circuit or dangerous malfunction

- Protect the handpiece against liquid penetration during cleaning and disinfection.

### Warnings

- The shield should be removed and cleaned/disinfected/sterilised. The handpiece, tips and charging base cannot be sterilised by autoclave. The device cannot tolerate high-level disinfection procedures.
- Intermediate-level disinfection is appropriate for the handpiece, tips and charging base.
- Do not clean or disinfect with chlorine bleach/sodium hypochlorite (corrosion of contacts) or Lysol® Brand I.C. Disinfectant Spray (cracking of charging base).
- Disconnect the power supply plug from the power outlet and charging base unit prior to cleaning/disinfection.

### Preparation before cleaning

Always remove the tip from the handpiece before processing. Use moderate force to pull the probe tip away from the handpiece.

### Cleaning

Do not use automated washer-disinfectors for reprocessing all the parts, as this will damage the components. They have to be manually cleaned.

#### *For handpiece, tip and charging base*

1. Remove the eye protection shield.
2. Remove the protective barrier sleeve and discard.
3. Scrub with an impregnated wipe or disposable towel soaked with an alcohol-based, tuberculocidal disinfectant, until it is free of visible residues.
4. Special note: use care when cleaning the mating surfaces of the probe tip and handpiece. Use only a damp, impregnated towel.
  - For the tip: vigorously scrub the area near the O-ring with a fresh wipe. Ensure that the liquid covers the O-ring and surrounding crevices. When cleaning the

mating surface, ensure that the cleaning agent only comes into contact with the sides that fit within the handpiece (with O-ring). Avoid applying cleaning agent to the electrical contacts on the bottom of the probe tip.

- For the handpiece mating cavity: use a fresh wipe to clean the mating groove directly below the surface. Use care to ensure that the cleaning agent is applied only to the top of the cavity interior. Ensure that only minimal cleaning agent enters the cavity that houses the electrical pins. Do not allow liquid to build up in the cavity around the contact pins. Immediately absorb any excess liquid with a dry disposable towel.
- For the battery and handpiece mating seam: use a fresh wipe to clean the mating groove. Remove all visible soiling, ensuring liquid penetrates all crevices. Use fresh wipes to rub liquid into the crevices. Do not allow any liquid to penetrate the casing. Discard used wipes. Additional wipes may be used.
- Do not remove the battery pack from the light handpiece. Do not attempt to disassemble the charging base.

5. Remove any cleaning solution residue with a damp cloth.

6. Allow the devices to air-dry for at least 5 min.

#### **For light shield**

1. Scrub with hot water and immerse the shield in a pH-neutral, phosphate-free cleaning detergent solution. Clean with a soft brush for at least 30 seconds until free of all visible contamination.
2. Rinse under running potable water.
3. Dry with a lint-free single-use cloth.

#### **Disinfection**

1. After cleaning, wipe all device surfaces with a new single-use cloth in combination with an alcohol-based, tuberculocidal, quaternary ammonium solution, with a contact time of 5 minutes, use according to the disinfectant solution manufacturer's instructions for use. Use a separate wipe for the tip and the handpiece. Ensure that the device and the disinfectant come into direct contact by pressing the wet wipes onto the device after half of the required contact time.
2. Use fresh wipes to disinfect the probe tip O-ring area, handpiece mating cavity, and battery/handpiece mating seam for the entire contact time. Carefully ensure cleaning agent is applied only to the top of the cavity interior. Ensure only minimal cleaning agent enters the cavity that houses the electrical pins. Immediately absorb any excess liquid with a dry disposable towel.
3. Wipe the device with a sterile, clean, lint-free cloth that is well dampened with deionised water for 30 seconds to remove all disinfectant. Pay special attention to all seams, especially around the probe tip/handpiece junction. Ensure that the cloth is damp with deionised water for the entire 30 seconds. Discard the used cloth and repeat rinsing with a new, second damp cloth for 30 seconds. Discard the second cloth and rinse once more with a new, third damp cloth for 30 seconds.

4. Wipe the device with a fourth dry, sterile, lint-free cloth to remove all liquid.
5. Allow the device to air-dry for at least 5 minutes.

### Packaging

No particular requirements.

### Sterilisation

#### ***For handpiece, tip and charging base***

- Must not be sterilised.
- Do not subject the components to steam autoclaving or liquid chemical sterilant solution immersion. This will cause damage to the components.

#### ***For light shield***

- After manual cleaning, and any optional disinfection or automated washer-disinfector cycle, steam autoclaving is required.
- Prevacuum steam sterilisation:
- Full cycle: 134 °C for 3 minutes 30 seconds.
- Follow the manufacturer's loading and operation cycle instructions.

### Drying

#### ***For handpiece, tip and charging base***

Wipe the devices dry with a sterile, clean, lint-free cloth. Allow the components to fully air-dry before storage.

#### ***For light shield***

Use the drying cycle of the autoclave, minimum 30 minutes. Allow the components to fully air-dry before storage.

### Maintenance, inspection and testing

- Carry out a visual inspection to ensure that all contamination has been removed. Visually inspect the power supply and cord for damage.
- Components that are damaged, worn or distorted, such as the O-rings, should be discarded and replaced.

### Storage

Store the components at room temperature, away from moisture or excessive humidity.

## 6. MAINTENANCE

### Light output monitoring

- Ensure that the LED aperture is clean and scratch-free; otherwise the light output will be reduced and may be insufficient for proper curing of the material.
- The light intensity should be checked frequently to ensure appropriate curing by using the radiometer incorporated into the charging base.

#### **Note:**

- The light radiometer can only be used to measure the curing light. For fluorescence detection light, the measured value is inaccurate.
- Before measuring, please turn the Rotary switch to Cure. Do not continue to use if the light output is below the intensity specified above.

### Battery

- When the battery light shows solid orange, the battery is charging. Once the battery is fully recharged, the battery light remains permanently green. The battery needs approximately 2 hours to fully recharge.
- If the battery pack needs to be replaced, simply remove the battery pack by pulling it from the main housing along its longitudinal axis.

### General maintenance

- A thin coating of petroleum jelly may be applied to the tip O-rings and charging base battery post as needed to facilitate insertion and removal.
- Inspect and replace worn or damaged O-rings as needed to maintain optimal performance.

## 7. TECHNICAL DATA, WARRANTY TERMS

|                                                   |                                                       |
|---------------------------------------------------|-------------------------------------------------------|
| AC supply                                         | AC100–240 V / 50–60 Hz                                |
| Charging base power input                         | SV, 1.5 A                                             |
| Operation ambient temperature                     | 5–40 °C                                               |
| Operation humidity                                | 20–80 %                                               |
| Transport and storage temperature                 | -10–55 °C                                             |
| Transport and storage humidity                    | < 93 % (non-condensing)                               |
| Battery performance                               | Recharge time: approximately 2 hours<br>3.7V, 300 mAh |
| Curing light intensity                            | More than 1,300 mW/cm <sup>2</sup>                    |
| Output peak wavelength range                      | 380–520 nm                                            |
| Handpiece dimensions<br>(including battery & tip) | Length 19.7 cm; Width 1.35 cm                         |
| Handpiece weight including battery pack           | 75.5 grammes                                          |

### Classifications

Type of protection against electric shock: Class II

Degree of protection against electric shock: Type B Applied Part

According to medical device directive: 1 (Rule 12) (IEC 60601)

### Device disposal

For environmental reasons, dispose of the device and battery according to local environmental guidelines or regulations.

### Warranty terms

- The manufacturer grants a 2-year warranty on all of the device's components except the battery. The battery is covered by a 1-year warranty. The warranty commences on the date of purchase. Within the warranty period, the manufacturer will eliminate free of charge any defects in the device resulting from faults in material or workmanship either by repairing or replacing parts or replacing the device in full at the manufacturer's discretion.

- Not covered by this warranty: damage caused by improper use (operation with incorrect current/voltage, unsuitable power point, breakage, cleaning via methods other than those recommended), normal wear and defects which have a negligible effect on the value or operation of the device.
- This warranty shall become void if repairs are undertaken by unauthorised persons.

## Symbol identification

|                                                                                  |                                                               |
|----------------------------------------------------------------------------------|---------------------------------------------------------------|
|  | Caution                                                       |
|  | Serial number                                                 |
|  | Type B applied part                                           |
|  | Keep dry                                                      |
|  | Vertical up                                                   |
|  | Direct current                                                |
|  | Special disposal of waste electrical and electronic equipment |

|                                                                                   |                             |
|-----------------------------------------------------------------------------------|-----------------------------|
|  | Warning                     |
|  | Do not reuse                |
|  | Fragile                     |
|  | Class II product            |
|  | Refer to instruction manual |
|  | Autoclave                   |

## Electromagnetic compatibility precaution

This information is required in accordance with the 4th edition of IEC 60601-1-2.

- The device requires special precautions regarding EMC and must be installed and put into service according to the EMC information provided in this manual.
- Portable and mobile RF communications equipment can affect the device.
- The use of accessories, transducers and cables other than those specified by manufacturer may result in increased emissions or decreased immunity of the device.
- The device should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
- In accordance with IEC 60601-1-2, no additional environmental operating conditions are required for normal use.

### Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

| Emissions test                                       | Compliance         | Electromagnetic environment – guidance                                                                                                                                                                                  |
|------------------------------------------------------|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions CISPR11                                 | Group 1            | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are unlikely to cause any interference in nearby electronic equipment.                                           |
| RF emissions CISPR11                                 | Class B            | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                     | Class A - Complies |                                                                                                                                                                                                                         |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Complies           |                                                                                                                                                                                                                         |

### Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

| Immunity test                               | IEC 60601 test level                                  | Compliance level                                      | Electromagnetic environment - guidance                                                                                                         |
|---------------------------------------------|-------------------------------------------------------|-------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact<br>± 2, ± 4, ± 8,<br>± 15<br>15 kV air | ± 8 kV contact<br>± 2, ± 4, ± 8,<br>± 15<br>15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |

|                                                                                     |                                                                                                                 |                                                                                                                 |                                                                                                                                                                                                                                                                              |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrical fast transient/burst<br>IEC 61000-4-4                                    | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for input/output lines                                          | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for input/output lines                                          | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                          |
| Surge<br>IEC 61000-4-5                                                              | $\pm 0.5$ , $+ 1$ kV line(s) to line(s)<br>$\pm 0.5$ , $\pm 1$ , $\pm 2$ kV line(s) to earth                    | $\pm 0.5$ , $+ 1$ kV line(s) to line(s)<br>$\pm 0.5$ , $\pm 1$ , $\pm 2$ kV line(s) to earth                    | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                          |
| Voltage dips, short interruptions and voltage variations on power<br>IEC 61000-4-11 | 0 % $U_T$ for 0.5 cycle<br>0 % $U_T$ for 1 cycle<br>70 % $U_T$ for 25/30 cycles<br>0 % $U_T$ for 250/300 cycles | 0 % $U_T$ for 0.5 cycle<br>0 % $U_T$ for 1 cycle<br>70 % $U_T$ for 25/30 cycles<br>0 % $U_T$ for 250/300 cycles | Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during mains power interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000 4-8                          | 30 A/m                                                                                                          | 30 A/m                                                                                                          | Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.                                                                                                                                                          |
| NOTE: $U_T$ is the a.c. mains voltage prior to application of the test level.       |                                                                                                                 |                                                                                                                 |                                                                                                                                                                                                                                                                              |

# **Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

| Immunity test                 | IEC 60601 test level                 | Compliance level   | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------------------------------|--------------------------------------|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000 4-6 | 3 V <sub>rms</sub> 150 kHz to 80 MHz | 3 V <sub>rms</sub> | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance<br><br>$d = ([3.5] \sqrt{p})/3$<br>$d = ([3.5] \sqrt{p})/380 \text{ MHz to } 800 \text{ MHz}$<br>$d = ([7] \sqrt{p})/3800 \text{ MHz to } 2,7 \text{ GHz}$<br><br>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |
| Radiated RF<br>IEC 61000 4-3  | 10 V/m<br>80 MHz to 2,7 GHz          | 10 V/m             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientation or relocating the [System Name here].
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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VOCO GmbH
**Lotus NL B.V.**
**Foshan COXO Medical Instrument Co.,Ltd**