Page 1 of 4

Revision: 1

## 1. Substance / Preparation and Company name

Product Name: riva self-cure hv capsules

Recommended use: For filling of cavitated teeth by dental professionals.

**Manufacturer / Supplier** 

SDI Limited SDI Inc.

5-9 Brunsdon Street, Bayswater 729 N.Route 83, Suite 315 Victoria, 3153, Australia Bensenville 60106 IL, USA

<u>Telephone</u>: <u>Telephone</u>:

+61 3 8727 7111 (Business hours) 630 238 8300 (Business hours)

Southern Dental Industries Ltd

Block 8, St Johns Court

Swords Road

SDI Brasil Indústria e Comércio Ltda
Rua Dr. Virgílio de Carvalho Pinto, 612
Pinheiros, São Paulo, 05415-020

Santry, Dublin 9, Ireland Brasil

<u>Telephone</u>: <u>Telephone</u>:

+353 1 886 9577 (Business Hours) +55 11 3031 1700 (Business Hours)

Emergency contact number: +61 3 8727 7111

## 2. Composition / Information on ingredients

Composition:	CAS No.	<u>Wt. %</u>
Compartment 1:		
Polyacrylic Acid Tartaric acid	9003-01-4 87-69-4	20-30 10-15
Compartment 2:		
Fluoro Aluminosilicate glass Polyacrylic Acid	- 9003-01-4	90-95 5-10

#### 3. Hazard Identification

Materials in the two compartments are irritants. However, under normal conditions of use, contact with these materials is generally not expected.

Risk Phrases: **36** Irritating to eyes.

**37** Irritating to respiratory system.

**38** Irritating to skin.

Safety Phrases: **25** Avoid contact with eyes.

In case of contact with eyes, rinse immediately with plenty of water and

seek medical advice.

28 After contact with skin, wash immediately with water and soap.

**MATERIAL SAFETY DATA SHEET** [according to GHS and NOHSC:2011(2003)]

**Product: RIVA SELF CURE HV CAPSULES** 

Date / Revised: 30.11.2010 Revision: 1

Page 1 of 4

#### 4. First Aid Measures

Eye (contact): Flush opened eye with running water for at least 5 minutes. Seek medical

attention.

Skin (contact): Remove contaminated clothing. Wash skin with plenty of water. If irritation

occurs seek medical attention.

Ingestion: Drink plenty of water. Seek medical attention if discomfort continues.

## 5. Fire Fighting Measures

Suitable extinguishing media: None required.

Unusual Fire and Explosion

Hazards: None.

Special protective equipment: No special protective equipment required.

#### 6. Accidental Release Measures

Personal precautions: Do not get into eyes, on skin or clothing.

Environmental precautions: Clear up spillages. Transfer to a container for disposal.

Methods for cleaning up: Wash the spillage area clean with water.

## 7. Handling and storage

Handling: No special measures required.

Storage: Store in a cool place at temperatures between 4°C and 20°C (39°F - 68°F).

### 8. Exposure controls and personal protection

Respiratory protection: None required.

Hand protection: Rubber, latex or PVC gloves.

Eye protection: Recommended.

General safety and hygiene measures: None.

**Product: RIVA SELF CURE HV CAPSULES** 

Date / Revised: 30.11.2010 **Revision: 1** 

Page 1 of 4

## Physical and chemical properties

Appearance: Smooth, pale-coloured paste

Odour: Slightly characteristic odour

Approx. 100°C (212°F) Boiling point:

Melting point: Not applicable

Specific gravity:  $2.0 \text{ g/cm}^3$ Flash point: Not applicable

Flammable: Not flammable Autoflammability: Not applicable

Explosive properties: Does not present an explosion hazard

Oxidizing properties: Not applicable Vapour pressure (@ 20°C): Approx. 23 hPa

Relative density: <1.0

Solubility: Miscible

## 10. Stability and Reactivity

Stable under normal conditions. Stability:

Conditions to avoid: Avoid excessive heat

Materials to avoid: None.

Hazardous decomposition products: None when used as directed.

Hazardous reactivity (polymerization): No hazardous reactions known.

### 11. Toxicological information

Contact with skin and eyes may cause irritation. Acute toxicity:

Eye (contact): Irritating to eyes

Skin (contact): Possible skin irritant

Inhalation: Not expected

Irritant not likely in small amounts Ingestion:

### 12. Ecological information

Self assessment: Not readily biodegradable. Avoid large quantities entering drains, surface water

or ground water.

**MATERIAL SAFETY DATA SHEET** [according to GHS and NOHSC:2011(2003)]

**Product: RIVA SELF CURE HV CAPSULES** 

Date / Revised: 30.11.2010 Revision: 1

Page 1 of 4

### 13. Disposal considerations

Dispose of in accordance with local official regulations.

# 14. Transport Information

Not classified as a Dangerous Good by road, rail, sea, or air.

### 15. Regulatory information

This product is regulated by

**TGA** 

Medical Devices Directive 93/42/EEC

**FDA** 

National regulations

#### 16. Other information

The information contained in the Material Safety Data Sheet is based on data considered to be accurate, however, no warranty is expressed or implied regarding the accuracy of the data or the results to be obtained from the use thereof.

+61 3 8727 7111

Prepared by: SDI Limited Phone Number:

3 -13 Brunsdon Street, Bayswater

Victoria, 3153, Australia

**Department issuing MSDS:** Research and Development

**Contact:** Operations Director