

020631

SAFETY DATA SHEET according to 91/155/EC

SDS n° 266 115 & 265 950

Date of issue: 06/07/2004

n.a. = not applicable - n.k. = not known



1 - IDENTIFICATION OF THE SUBSTANCE / PREPARATION OF THE COMPANY / UNDERTAKING

Product details:

Trade name: SALVIZOL E.D.T.A.

Manufacturer / Supplier: PRODUITS DENTAIRE PIERRE ROLLAND
ZI PHARE / BP 216
33708 MERIGNAC CEDEX FRANCE

Telephone: + 33 (0) 556 34 06 07
Fax: + 33 (0) 556 34 92 92

Emergency Information: Call + 33 (0) 556 34 06 07

2 - COMPOSITION / INFORMATION ON INGREDIENTS

Chemical characterization:

Tetracetic diamine ethylene acid (CAS NR 60-00-4), concentration: 8 per cent, Xi, R36

Benzalkonium chloride (CAS NR 63449-41-2), concentration: 0,3 per cent, C, R21/22, R34

3 - HAZARDS IDENTIFICATION

No health or environment hazards when the material is handled in accordance with official regulations.

4 - FIRST-AID MEASURES

General information:

After inhalation: Breathe fresh air

After skin contact: Wash immediately with water and soap

After eye contact: Rinse immediately and abundantly for 10 minutes keeping the eyelids opened. Seek medical advice.

After swallowing: Seek medical advice immediately.

5 - FIRE-FIGHTING MEASURES

Suitable extinguishing media: All fire-extinguishers are suitable.

Extinguishing media which must not be used for safety reasons: None according to our knowledge.

6 - ACCIDENTAL RELEASE MEASURES

Environmental precautions: Do not discard into sewers or streams. Do not discard with household waste.

Methods for cleaning up: Pick up with appropriate products and store for a late disposal in accordance with legislation in force.

7 - HANDLING AND STORAGE

Handling: Keep out of the reach of children. Do not swallow. Avoid projections on the skin and in the eyes. Do not breathe vapours. Well close after use.

Storage: Store in its original packaging at a temperature not exceeding 25 °C and protected from heat. Keep the jar well closed between each use.

8 - EXPOSURE CONTROLS / PERSONAL PROTECTION

Hand and skin protection: Wear appropriate gloves, goggles and clothes.

Hygiene measures: Do not eat, drink or smoke when using.

